Question 1

A. 100ml of 0.0250 M solution of codeine (Mw = 299.37 g/mol)

```
M=n/V so, n=M*V=0.0250 M*0.1L=0.0025 moles of codeine.
Amount= n*Mw=(0.0025 moles*299.37g)/1 mol= 0.7484g of codeine.
```

B. C1V1=C2V2

C1= 0.0250 M

V1= x

C2 = 0.001M

V2= 0.250 L

V1=(C2V2)/C1= (0.001*0.250)/ 0.0250= **0.01L= 10mL**

Question 2

- A) 100 ml of 1 mg/mL solution= 100*1= 100 mg
- B) 10 ml of 10 mg/mL =10*10= 100 mg
- C) 20 ml of 20 mg/L solution = 0.02L*20 mg/L = 0.4 mg
- D) 50 ml of a 10 mg/mL solution = 50*10=500mg

D contains the highest amount of drug

Question 3

```
A = \epsilonlc so c= A/(\epsilonl)
A=0.345
\epsilon= 569 L/(mol*cm)
l= 2cm
```

c (mol*L-1)= 0.345/ (569 L/(mol*cm) *2cm) = $3.03*10^{-4}$ mol*L⁻¹ As 4 liters = $3.03*10^{-4}$ mol*L⁻¹ *4L= 0.001213 mol

Amount= c*Mw= 0.001213 mol *185g/mol = **0.224** g

Question 4

 $Na_2CO_3 = 250 \text{ mL*}0.005 \text{mmol/mL} = 1.25 \text{ mmol/}2 = 0.625 \text{ mmol} = 6.25 \text{x} 10^{-4} \text{ mol *} 105.99 \text{ g} = 0.066 \text{g} Na_2CO_3$

Dissolving 0.066g Na₂CO₃ in water and diluting 250 mL

Question 5

HCL + NaOH -> NaCL + H20

0.105M x 0.005 L 0.0186 L

M=n/V so, n= M*V= 0.105*0.005= 5.25x10-4 mol HCL NaOH 1:1 HCl so, 5.25x10-4 mol HCL= 5.25x10-4 mol NaOH M NaOH = 5.25x10-4 mol NaOH/0.0186= **0.028M NaOH**

Question 6

Kd= 4.0 Vo= 50 mL Va= 100mL

$$F = \frac{Kd}{Kd + \left(\frac{Va}{Vo}\right)} = \frac{4}{4 + \left(\frac{100}{50}\right)} = \frac{4}{6} = 0.6667 = 66.67\%$$

Question 7

1 mg in 100 ml iron 40% light 0.5 mg iron in 100 ml

A=-log(I/Io)=-log(40/100)= -log(0.4)= 0.398 0.398*0.5= 0.199 T= 10 -0.199 = 0.6324 = **63.24**%

Question 8

Accuracy defines how close the result is to the true value. Precision is a measure of repeatability, quoted as absolute standard deviation (S) or relative standard deviation (RSD).

Average= 3.9875g

Theoretical value= 5mL of 0.78 g/mL solution= 3.9g

$$\text{Relative error=} \left(\frac{experimental\ value-\ theoretical\ value}{theoretical\ value} \right) *\ 100 = \left(\frac{3.9875g - 3.9g}{3.9875g} \right) *\ 100 = 2.24\%$$

The lower the relative error is, the more accurate the pipette is.

$$S=[\sum (X i-X)^2/(N-1)]^{1/2}$$

$$\sum (X \text{ i}-X)^2 = 2(3.9-3.9875)^2 + 3(4-3.9875)^2 + (3.8-3.9875)^2 + (4.2-3.9875)^2 + (4.1-3.9875)^2 = 0.10875$$

 $S = \left[\sum (X \text{ i}-X)^2/(N-1)\right]^{1/2} = (0.10875/7)^{\frac{1}{2}} = 0.125$

RSD= (s/average)*100 = (0.125/3.9875)*100= 3.135%

Looking at the low value of the results, the pipettes are acceptable regarding their accuracy and precision.

Question 9

```
H_2SO_4 + 2NaOH -> Na_2SO_4 + 2H_2O

2OmL 14.5 mL

x 0.150M
```

M=n/V so, n= M*V=0.150*0.0145= 0.002175 moles NaOH H_2SO_41 : 2 NaOH -> 0.002175 mole

 $0.002175 \text{ moles NaOH} -> 1 \text{mol H}_2 \text{SO}_4 / 2 \text{mol NaOH} = 0.002175 \text{ moles NaOH}/2 = 1.0875 \times 10^{-3} \text{ moles H}_2 \text{SO}_4$

 $M = 1.0875 \times 10^{-3} \text{ moles } H_2SO_4/0.0206 \text{ L} = 5.279 \times 10^{-2} \text{M } H_2SO_4$

Question 10

- i. Concentration (in mg/L) of this standard iron solution c= 0.1483 g in 250 mL= 0.5932 g/L =**593.2 mg/L**
- ii. 1 mL of 0.5 mg/L stock solution

 $C_1 = 0.5 \text{ mg/L}$

 $V_1=1 \text{ mL} = 0.001$

 $C_2 = x$

 $V_2 = 100 \text{ mL} = 0.1 \text{L}$

$$C1 * V1 = C2 * V2 \rightarrow C2 = \frac{C1 * V1}{V2} = \frac{0.001 * 0.5}{0.1} = \mathbf{0.005} \, mg/L$$

3 mL of 0.5 mg/L stock solution

 $C_1 = 0.5 \text{ mg/L}$

 $V_1 = 3 \text{ mL} = 0.003 \text{ L}$

 $C_2 = x$

 $V_2 = 100 \text{ mL} = 0.1 \text{L}$

$$C1 * V1 = C2 * V2 \rightarrow C2 = \frac{C1 * V1}{V2} = \frac{0.003 * 0.5}{0.1} = \mathbf{0.015} \ mg/L$$

5 mL of 0.5 mg/L stock solution

 $C_1 = 0.5 \text{ mg/L}$

 $V_1 = 5 \text{ mL} = 0.005 \text{ L}$

 $C_2 = x$

 $V_2 = 100 \text{ mL} = 0.1 \text{L}$

$$C1 * V1 = C2 * V2 \rightarrow C2 = \frac{C1 * V1}{V2} = \frac{0.005 * 0.5}{0.1} = 0.025 \, mg/L$$

7 mL of 0.5 mg/L stock solution

 $C_1 = 0.5 \text{ mg/L}$

 $V_1 = 7 \text{ mL} = 0.007 \text{ L}$

 $C_2 = x$

 $V_2 = 100 \text{ mL} = 0.1 \text{L}$

$$C1 * V1 = C2 * V2 \rightarrow C2 = \frac{C1 * V1}{V2} = \frac{0.007 * 0.5}{0.1} = \mathbf{0.035} \ mg/L$$

10 mL of 0.5 mg/L stock solution


 $C_1 = 0.5 \text{ mg/L}$

 $V_1 = 10 \text{ mL} = 0.010 \text{ L}$

 $C_2 = x$

 $V_2 = 100 \text{ mL} = 0.1 \text{L}$

$$C1 * V1 = C2 * V2 \rightarrow C2 = \frac{C1 * V1}{V2} = \frac{0.010 * 0.5}{0.1} = \mathbf{0.05} \ mg/L$$



iii.

Line equation: y=0.1738 L*mg-1*Abs-1x-0.029 Abs

iv. Concentration of the Unknown sample= (0.375 Abs + 0.029 Abs)/ 0.1738 L*mg-1*Abs-1 =**2.325 mg/L**

Y=mx+c

y-c=mx

X=(y-c)/m

Question 11

$$Kd = 5.32$$

$$Va = 100 \text{ ml}$$

$$Vo = 30 \text{ ml}$$

$$N = 3$$

$$F = 1 - \left(\frac{Va}{Kd * Vo + Va}\right)^3 = 1 - \left(\frac{100}{5.32 * 30 + 100}\right)^3 = 0.9428 = 94.28\%$$

Question 12



Question 13

Solution 1

 $C_1 = 20 \mu g/L$

 $V_1=x$

 $C_2 = 10 \mu g/L$

 $V_2 = 20 \text{ mL}$

$$C1 * V1 = C2 * V2 \rightarrow V1 = \frac{C2 * V2}{C1} = \frac{10 * 20}{20} = 10mL$$

Solution 2

 $C_1 = 20 \mu g/L$

 $V_1=x$

 $C_2 = 5 \mu g/L$

V₂= 20 mL

$$C1 * V1 = C2 * V2 \rightarrow V1 = \frac{C2 * V2}{C1} = \frac{20 * 5}{20} = 5 \ mL$$

Solution 3

 $C_1 = 20 \mu g/L$

 $V_1=x$


 $C_2 = 2.5 \mu g/L$

 $V_2 = 20 \text{ mL}$

$$C1 * V1 = C2 * V2 \rightarrow V1 = \frac{C2 * V2}{C1} = \frac{20 * 2.5}{20} = 2.5 mL$$

We should take the appropriate volume of the stock solution with the right pipette. A 5mL pipette should be used for solutions 1 and 2. For solution 3, a 2mL pipette would be more accurate.

Question 14



Line equation: y=0.2485x-0.2092

Dilution factor= 250

Amount of drug in tablet (0.386 Abs+0.2092 Abs)/0.2485 L*mg-1*Abs-1= 2.395 mg/L*250= 598.75 mg/L

There is 0.599g of drug in the tablet (1.355g), so 0.599g/1.355g= 0.4419g*100= 44.19%

There is **44.19%** of drug in the tablet.

Question 15

Va= 5mL Vo= 10 mL pH= 2 Ka= 3.24×10^{-4} Kd= 5

$$[H^+]$$
= pH 2= 10^{-2} = $1x10^{-2}$

$$D = \frac{Kd*[H+]}{[H+]+Ka} = \frac{5*1x10-2}{1x10-2+3.24 \times 10-4} = 4.84$$

$$F = \frac{Kd}{Kd + (\frac{Va}{Vo})} = \frac{4.84}{4.84 + (\frac{5}{10})} = 0.9064 = \mathbf{90.64}\%$$